



LET'S WORK TOGETHER

Here to help with
your essentials

LED LIGHT SOURCE

LED - 1500

LED Light Source designed to provide
light for endoscopic procedures.



Discover more of our product at: www.divinemeditech.in



LED LIGHT SOURCE

LED - 1500

LED Light Source designed to provide light for endoscopic procedures.



Salient Features

- Focus coherent lens technique for uniform brightness
- Special heat filter
- Continuous variable light intensity 1-100 %
- Extremely high optical output of light
- LCD Display of light intensity and lamp life
- Low noise level enabling quiet operation
- Optional remote-control facility for essential functions
- Standby mode
- Compact and Light weight construction

Technical Specification

- | | |
|-------------------------|--|
| • Type of LED | : Large area LED chip, ROHS, REACH compliant for light |
| • Illumination | : Equivalent to 300w xenon light |
| • Color Temperature | : >6500K |
| • Standard Connection | : Storz type, other types on request |
| • Classification | : Class 1 / Class A Medical Device |
| • Ambient Temperature | : +10° C to +40° C during storage |
| • Dimension (W x H x D) | : 35 x 13 x 37 cm |
| • Power Supply | : 100 / 240 + 10% VAC, 50/60 Hz |
| • Fuses | : 2 x 3 Amp, 250 VAC |
| • Atmospheric Pressure | : 86 kPa-106kPa |